

C1Q-Polyclonal-Antibody-HRP-Conjugated

Catalog No: #C90678H



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	C1Q-Polyclonal-Antibody-HRP-Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC
Species Reactivity	Hu Ms Rat
Conjugates	HRP
Target Name	C1Q
Other Names	C1Q
Accession No.	Uniprot:P02745Gene ID:712
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:100-500 IHC-F=1:50-200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.