

Anti-VKORC1 Antibody

Catalog No: #62003



Package Size: #62003-1 50ul #62003-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Anti-VKORC1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Human
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the N-terminal region of human VKORC1. The exact sequence is proprietary.
Calculated MW	Predicted band size: 18 kD; Observed band size: 15 kD
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

WB 1αO 500 - 1αO 1000

Images



Western blot analysis of VKORC1 expression in human lung (A) whole cell lysates

Note: This product is for in vitro research use only and is not intended for use in humans or animals.