

VEGFR2 Polyclonal Antibody

Catalog No: #C90236



Package Size: #C90236-1 50ul #C90236-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	VEGFR2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC-P, Flow-Cyt
Species Reactivity	Human, Mouse, Rat
Calculated MW	147
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P □ 1:100-500

Flow-Cyt □ 1 μg/Test

Note: This product is for in vitro research use only and is not intended for use in humans or animals.