

Human Complement C1q Tumor Necrosis Factor Related protein 4 (C1QTNF4) ELISA Kit

Catalog No: #EK20761



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #EK20761-1 48T #EK20761-2 96T

Description

Product Name	Human Complement C1q Tumor Necrosis Factor Related protein 4 (C1QTNF4) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.156-10 ng/mL
Sensitivity:0.067ng/mL
Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.