

B-Myb Polyclonal Antibody

Catalog No: #62006



Package Size: #62006-1 50ul #62006-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	B-Myb Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IP IF ELISA
Species Reactivity	Human Mouse Rat Rice
Specificity	B-Myb Polyclonal Antibody detects endogenous levels of B-Myb protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human B-Myb.
Target Name	MYBL2
Other Names	MYBL2;BMYB;Myb-related protein B;B-Myb;Myb-like protein 2
Calculated MW	80kD
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

WB 1:500 - 1:2000


IF 1:200 - 1:1000.

ELISA: 1:20000

Images



Western blot analysis of lysates from A549 cells, using B-Myb Antibody. The lane on the right is blocked with the synthesized peptide.



Immunofluorescence analysis of HeLa cells, using B-Myb Antibody. The picture on the right is blocked with the synthesized peptide.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.