

APC anti-human CD273

Catalog No: #FC0104



Package Size: #FC0104-1 25Tests #FC0104-2 100Tests

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	APC anti-human CD273
Brief Description	mouse monoclonal antibody
Host Species	mouse
Clonality	Monoclonal
Clone No.	s-24F.10C12
Isotype	mouse IgG2a, κ
Applications	FC
Species Reactivity	human
Conjugates	APC
Target Name	CD273
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



Human monocyte-derived dendritic cells stained with APC anti-human CD273 and PE anti-human CD11c

Note: This product is for in vitro research use only and is not intended for use in humans or animals.