

Biotin anti-human CD326

Catalog No: #FC0145



Package Size: #FC0145-1 25ug #FC0145-2 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	Biotin anti-human CD326
Brief Description	mouse monoclonal antibody
Host Species	mouse
Clonality	Monoclonal
Clone No.	s-VU-1D9
Isotype	mouse IgG1
Applications	FC
Species Reactivity	human
Conjugates	Biotin
Target Name	CD326
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



Two mixed HT-29 cell line stained with Biotin anti-human CD326 (colored histogram) or Biotin mouse IgG1 (gray histogram), followed by SAV-PE.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.