

APC-Cy7 anti-human CD25

Catalog No: #FC0449



Package Size: #FC0449-1 25Tests #FC0449-2 100Tests

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	APC-Cy7 anti-human CD25
Brief Description	mouse monoclonal antibody
Host Species	mouse
Clonality	Monoclonal
Clone No.	s-MEM-181
Isotype	mouse IgG1
Applications	FC
Species Reactivity	human
Conjugates	APC-Cy7
Target Name	CD25
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



PHA-stimulated (3 days) human peripheral blood lymphocytes stained with APC-Cy7 anti-human CD25 (colored histogram) or APC-Cy7 mouse IgG1 (gray histogram)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.