

PE anti-human CD115

Catalog No: #FC0582



Package Size: #FC0582-1 25Tests #FC0582-2 100Tests

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	PE anti-human CD115
Brief Description	rat monoclonal antibody
Host Species	rat
Clonality	Monoclonal
Clone No.	s-2-s-5-4
Isotype	rat IgG1,k
Applications	FC
Species Reactivity	human
Conjugates	PE
Target Name	CD115
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



Human peripheral blood monocytesstained with PE anti-human 115 (colored histogram) or PE Rat IgG1 κ(gray histogram)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.