

PE anti-human/mouse/rat CD278 (ICOS)

Catalog No: #FC0610



Package Size: #FC0610-1 25Tests #FC0610-2 100Tests

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


| | |
|--------------------|---|
| Product Name | PE anti-human/mouse/rat CD278 (ICOS) |
| Brief Description | Armenian Hamster Monoclonal Antibody |
| Host Species | Armenian Hamster |
| Clonality | Monoclonal |
| Clone No. | s-C398.s- |
| Isotype | Armenian Hamster IgG |
| Applications | FC |
| Species Reactivity | human/mouse/rat |
| Conjugates | PE |
| Target Name | CD278 (ICOS) |
| Formulation | 0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4 |
| Storage | Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL. |

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



PHA-stimulated (3 days) PBMC stained with PE anti-human CD278 (colored histogram) or PE Armenian hamster IgG (gray histogram)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.