

PerCP-Cy5.5 anti-human/mouse CD11b

Catalog No: #FC0680



Package Size: #FC0680-1 25Tests #FC0680-2 100Tests

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


| | |
|--------------------|---|
| Product Name | PerCP-Cy5.5 anti-human/mouse CD11b |
| Brief Description | rat monoclonal antibody |
| Host Species | rat |
| Clonality | Monoclonal |
| Clone No. | s-M1/70 |
| Isotype | rat IgG2b |
| Applications | FC |
| Species Reactivity | human/mouse |
| Conjugates | PerCP-Cy5.5 |
| Target Name | CD11b |
| Formulation | 0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4 |
| Storage | Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL. |

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



Human peripheral blood granulocytes stained with PerCP-Cy5.5 anti-human /mouse CD11b (colored histogram) or PerCP-Cy5.5rat IgG2b (gray histogram)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.