

# PerCP-Cy5.5 anti-human CD127

Catalog No: #FC0699



Package Size: #FC0699-1 25Tests #FC0699-2 100Tests

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	PerCP-Cy5.5 anti-human CD127
Brief Description	mouse monoclonal antibody
Host Species	mouse
Clonality	Monoclonal
Clone No.	s-A019D5
Isotype	mouse IgG1,k
Applications	FC
Species Reactivity	human
Conjugates	PerCP-Cy5.5
Target Name	CD127
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

## Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

## Images



Human peripheral blood CD4+ T cells stained with PerCP-Cy5.5 anti-human CD127 and PE anti-human CD25

Note: This product is for in vitro research use only and is not intended for use in humans or animals.