

Pacific Blue anti-human/mouse/rat CD278 (ICOS)

Catalog No: #FC0742



Package Size: #FC0742-1 25Tests #FC0742-2 100Tests

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	Pacific Blue anti-human/mouse/rat CD278 (ICOS)
Brief Description	Armenian Hamster Monoclonal Antibody
Host Species	Armenian Hamster
Clonality	Monoclonal
Clone No.	s-C398.s-
Isotype	Armenian Hamster IgG
Applications	FC
Species Reactivity	human/mouse/rat
Conjugates	Pacific Blue
Target Name	CD278 (ICOS)
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



PHA-stimulated (3 days) PBMC stained with Pacific Blue anti-human CD278 (colored histogram) or Pacific Blue Armenian hamster IgG (gray histogram)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.