

# PE anti-human FOXP3

Catalog No: #FC0801

Package Size: #FC0801-1 25Tests #FC0801-2 100Tests

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	PE anti-human FOXP3
Brief Description	rat monoclonal antibody
Host Species	rat
Clonality	Monoclonal
Clone No.	s-PCH101
Isotype	rat IgG2a, k
Applications	FC
Species Reactivity	human
Conjugates	PE
Target Name	FOXP3
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

## Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

## Images



Human peripheral blood lymphocytes was fixed permeabilized and intracellular stained with PE anti-human FOXP3 and APC anti-human CD25

Note: This product is for in vitro research use only and is not intended for use in humans or animals.