

# APC anti-mouse FOXP3

Catalog No: #FC0914



Package Size: #FC0914-1 25Tests #FC0914-2 100Tests

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	APC anti-mouse FOXP3
Brief Description	rat monoclonal antibody
Host Species	rat
Clonality	Monoclonal
Clone No.	s-MF-14
Isotype	rat IgG2b, κ
Applications	FC
Species Reactivity	mouse
Conjugates	APC
Target Name	FOXP3
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.


## Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

## Images





---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.