

PE-Cy5 anti-mouse CD3ε

Catalog No: #FC0929



Package Size: #FC0929-1 25Tests #FC0929-2 100Tests

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	PE-Cy5 anti-mouse CD3ε
Brief Description	Armenian Hamster Monoclonal Antibody
Host Species	Armenian Hamster
Clonality	Monoclonal
Clone No.	s-145-2C11
Isotype	Armenian Hamster IgG
Applications	FC
Species Reactivity	mouse
Conjugates	PE-Cy5
Target Name	CD3ε
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.