

# PE-Cy5 anti-mouse CD8 $\alpha$

Catalog No: #FC0932



Package Size: #FC0932-1 25Tests #FC0932-2 100Tests

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | PE-Cy5 anti-mouse CD8 $\alpha$                                    |
| Brief Description  | rat monoclonal antibody   |
| Host Species       | rat   |
| Clonality          | Monoclonal  |
| Clone No.          | s-53-6.7  |
| Isotype            | rat IgG2a, $\kappa$   |
| Applications       | FC  |
| Species Reactivity | mouse   |
| Conjugates         | PE-Cy5  |
| Target Name        | CD8 $\alpha$  |
| Formulation        | 0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4 |
| Storage            | Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.            |

## Application Details

Test: 5  $\mu$ l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.