

# FITC anti-mouse CD30

Catalog No: #FC0981



Package Size: #FC0981-1 25Tests #FC0981-2 100Tests

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
 Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	FITC anti-mouse CD30
Brief Description	Armenian Hamster Monoclonal Antibody
Host Species	Armenian Hamster
Clonality	Monoclonal
Clone No.	s-mCD30.1
Isotype	Armenian Hamster IgG1
Applications	FC
Species Reactivity	mouse
Conjugates	FITC
Target Name	CD30
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

## Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

## Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.