

FITC anti-mouse H-2Kd (MHC Class I)

Catalog No: #FC1034

Package Size: #FC1034-1 25Tests #FC1034-2 100Tests

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	FITC anti-mouse H-2Kd (MHC Class I)
Brief Description	mouse monoclonal antibody
Host Species	mouse
Clonality	Monoclonal
Clone No.	s-A4C8.1-Do9
Isotype	mouse IgG1,k
Applications	FC
Species Reactivity	mouse
Conjugates	FITC
Target Name	H-2Kd (MHC Class I)
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.