

PE anti-mouse TCR $\gamma\delta$

Catalog No: #FC1188



Package Size: #FC1188-1 25Tests #FC1188-2 100Tests

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	PE anti-mouse TCR $\gamma\delta$
Brief Description	Armenian Hamster Monoclonal Antibody
Host Species	Armenian Hamster
Clonality	Monoclonal
Clone No.	s-UC7-13D5
Isotype	Armenian Hamster IgG
Applications	FC
Species Reactivity	mouse
Conjugates	PE
Target Name	TCR $\gamma\delta$
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 5 μ l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.