

# Pacific Blue anti-mouse Ly-6C

Catalog No: #FC1216



Package Size: #FC1216-1 25Tests #FC1216-2 100Tests

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Pacific Blue anti-mouse Ly-6C                                     |
| Brief Description  | rat monoclonal antibody   |
| Host Species       | rat   |
| Clonality          | Monoclonal  |
| Clone No.          | s-HK1.4   |
| Isotype            | rat IgG2c, κ  |
| Applications       | FC  |
| Species Reactivity | mouse   |
| Conjugates         | Pacific Blue  |
| Target Name        | Ly-6C   |
| Formulation        | 0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4 |
| Storage            | Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.            |

## Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.