

# PE Mouse IgG2a $\kappa$ isotype Control

Catalog No: #FC1258



Package Size: #FC1258-1 25Tests #FC1258-2 100Tests

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PE Mouse IgG2a $\kappa$ isotype Control
Clonality	Monoclonal
Clone No.	s-MOPC-173
Applications	FC
Species Reactivity	Mouse
Conjugates	PE
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

## Application Details

Test: 5  $\mu$ l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically,  $1 \times 10^6$  cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.