

# Alexa Fluor 488 Armenian Hamster IgG isotype Control

Catalog No: #FC1312



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #FC1312-1 25Tests #FC1312-2 100Tests

## Description

Product Name	Alexa Fluor 488 Armenian Hamster IgG isotype Control
Clonality	Monoclonal
Clone No.	s-PIP
Applications	FC
Species Reactivity	Armenian Hamster
Conjugates	Alexa Fluor 488
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

## Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.