

Alexa Fluor 488 Mouse IgG1 κ isotype Control

Catalog No: #FC1319



Package Size: #FC1319-1 25Tests #FC1319-2 100Tests

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Alexa Fluor 488 Mouse IgG1 κ isotype Control
Clonality	Monoclonal
Clone No.	s-IB7.11
Applications	FC
Species Reactivity	Mouse
Conjugates	Alexa Fluor 488
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 5 μ l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1×10^6 cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.