

Alexa Fluor 647 Mouse IgM isotype Control

Catalog No: #FC1324



Package Size: #FC1324-1 25Tests #FC1324-2 100Tests

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Alexa Fluor 647 Mouse IgM isotype Control |
| Clonality | Monoclonal |
| Clone No. | s-7C2C6C12 |
| Applications | FC |
| Species Reactivity | Mouse |
| Conjugates | Alexa Fluor 647 |
| Formulation | 0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4 |
| Storage | Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL. |

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.