

# PE-Cy7 Rat IgG2a isotype Control

Catalog No: #FC1345



Package Size: #FC1345-1 25Tests #FC1345-2 100Tests

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | PE-Cy7 Rat IgG2a isotype Control                                  |
| Clonality          | Monoclonal  |
| Clone No.          | s-1-1   |
| Applications       | FC  |
| Species Reactivity | Rat   |
| Conjugates         | PE-Cy7  |
| Formulation        | 0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4 |
| Storage            | Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.            |

## Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.