

FITC anti-human/mouse CD44

Catalog No: #FC1370



Package Size: #FC1370-1 25Tests #FC1370-2 100Tests

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	FITC anti-human/mouse CD44
Brief Description	rat monoclonal antibody
Host Species	rat
Clonality	Monoclonal
Clone No.	s-IM7
Isotype	rat IgG2b, κ
Applications	FC
Species Reactivity	human/mouse
Conjugates	FITC
Target Name	CD44
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 5 μ l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1×10^6 cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



Human peripheral blood lymphocytes stained with FITC anti-human CD44 (colored histogram) or FITC rat IgG2b kgray histogram)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.