

Alexa Fluor 488 anti-mouse Foxp3

Catalog No: #FC1385



Package Size: #FC1385-1 25Tests #FC1385-2 100Tests

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	Alexa Fluor 488 anti-mouse Foxp3
Brief Description	rat monoclonal antibody
Host Species	rat
Clonality	Monoclonal
Clone No.	s-SB168a
Isotype	rat IgG2a
Applications	FC
Species Reactivity	mouse
Conjugates	Alexa Fluor 488
Target Name	Foxp3
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.


Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Images





Note: This product is for in vitro research use only and is not intended for use in humans or animals.