

Purified anti-human CD8

Catalog No: #FC1433



Package Size: #FC1433-1 100ug #FC1433-2 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|---------------------------|
| Product Name | Purified anti-human CD8 |
| Brief Description | mouse monoclonal antibody |
| Host Species | mouse |
| Clone No. | s-UCHT-4 |
| Species Reactivity | human |
| Storage | Store at 4°C. |

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.