

# Purified anti-human CD11b

Catalog No: #FC1440

Package Size: #FC1440-1 100ug #FC1440-2 1mg



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Purified anti-human CD11b
Brief Description	mouse monoclonal antibody
Host Species	mouse
Clone No.	s-LM2/1.6.11
Species Reactivity	human
Storage	Store at 4°C.

## Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.