

Purified anti-human CD13

Catalog No: #FC1446



Package Size: #FC1446-1 100ug #FC1446-2 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Purified anti-human CD13
Brief Description	mouse monoclonal antibody
Host Species	mouse
Clone No.	s-46A11E3
Species Reactivity	human
Storage	Store at 4°C.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.