

Purified anti-human IFN-r

Catalog No: #FC1557



Package Size: #FC1557-1 100ug #FC1557-2 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Purified anti-human IFN-r
Brief Description	mouse monoclonal antibody
Host Species	mouse
Clone No.	s-GAMMA3-11.1
Species Reactivity	human
Storage	Store at 4°C.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.