

Purified anti-human/mouse integrin β 7

Catalog No: #FC1565



Package Size: #FC1565-1 100ug #FC1565-2 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Purified anti-human/mouse integrin β 7
Brief Description	rat monoclonal antibody
Host Species	rat
Clone No.	s-FIB21
Species Reactivity	human/mouse
Storage	Store at 4°C.

Application Details

Test: 5 μ l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1×10^6 cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.