

# Purified anti-mouse CD16.2

Catalog No: #FC1589

Package Size: #FC1589-1 100ug #FC1589-2 1mg



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Purified anti-mouse CD16.2
Brief Description	Armenian Hamster Monoclonal Antibody
Host Species	Armenian Hamster
Clone No.	s-9E9
Species Reactivity	mouse
Storage	Store at 4°C.

## Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.