

# Purified anti-mouse CD19

Catalog No: #FC1590

Package Size: #FC1590-1 100ug #FC1590-2 1mg



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                          |
|--------------------|--------------------------|
| Product Name       | Purified anti-mouse CD19 |
| Brief Description  | rat monoclonal antibody  |
| Host Species       | rat                      |
| Clone No.          | s-1D3                    |
| Species Reactivity | mouse                    |
| Storage            | Store at 4°C.            |

## Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.