

Porcine Angiotensin 1-7 (Ang1-7) ELISA Kit

Catalog No: #EK20753



Package Size: #EK20753-1 #EK20753-2

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Porcine Angiotensin 1-7 (Ang1-7) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Porcine

Application Details

Detect Range:78.1-5000pg/mL

Sensitivity:29pg/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.