

# Porcine Angiotensin II (AngII) ELISA Kit

Catalog No: #EK20757



Package Size: #EK20757-1 #EK20757-2

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Porcine Angiotensin II (AngII) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Porcine

## Application Details

Detect Range:7.81-500 pg/mL

Sensitivity:2.3pg/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.