

Mouse N-Acylsphingosine Amidohydrolase 1 (ASAH1) ELISA Kit

Catalog No: #EK20765

Package Size: #EK20765-1 48T #EK20765-2 96T



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse N-Acylsphingosine Amidohydrolase 1 (ASAH1) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse
Target Sequence	Detection Method:sandwich

Application Details

Detect Range:0.312-20ng/mL

Sensitivity:0.122ng/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.