

Human Erythropoietin (EPO) ELISA Kit

Catalog No: #EK20772



Package Size: #EK20772-1 48T #EK20772-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human Erythropoietin (EPO) ELISA Kit
Applications	ELISA
Species Reactivity	Human
Target Sequence	Detection Method:sandwich

Application Details

Detect Range:31.2-2000pg/mL

Sensitivity:12.7pg/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.