

Vinculin mouse mAb

Catalog No: #62029



Package Size: #62029-1 50ul #62029-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vinculin mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Purification	Protein G
Applications	WBoO IFoO ELISA
Species Reactivity	HumanO MouseO Rat
Specificity	This antibody detects endogenous levels of Vinculin protein.
Immunogen Description	Synthesized peptide derived from human Vinculin.
Target Name	VCL
Other Names	VCL;Vinculin;Metavinculin
Calculated MW	124kD
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000.

IF 1:100-500.

ELISA 1:1000-5000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.