

PINK1 (Phospho-S228) Antibody

Catalog No: #62030



Package Size: #62030-1 50ul #62030-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	PINK1 (Phospho-S228) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Immunogen Description	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S228 of human PINK1 protein.
Other Names	Serine/threonine-protein kinase PINK1 mitochondrial; BRPK; PTEN-induced putative kinase protein 1
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	store at -20°C for one year

Application Details

WB 1:500 - 1000

IHC 1:50 - 100

Note: This product is for in vitro research use only and is not intended for use in humans or animals.