

Mouse Squamous Cell Carcinoma Antigen (SCCA) ELISA Kit

Catalog No: #EK20788

Package Size: #EK20788-1 48T #EK20788-2 96T



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Squamous Cell Carcinoma Antigen (SCCA) ELISA Kit
Target Sequence	Detection Method:sandwich

Application Details

Detect Range:62.5-4000pg/mL

Sensitivity:23.6pg/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.