

Porcine Beta-2-Microglobulin (b2M) ELISA Kit

Catalog No: #EK20789



Package Size: #EK20789-1 48T #EK20789-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name Porcine Beta-2-Microglobulin (b2M) ELISA Kit

Target Sequence Detection Method:sandwich

Application Details

Detect Range:3.12-200ng/mL

Sensitivity:1.27ng/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.