

UBAP2 Polyclonal Antibody

Catalog No: #62032

Package Size: #62032-1 50ul #62032-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | UBAP2 Polyclonal Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Applications | WB,ELISA |
| Species Reactivity | Human,Mouse |
| Specificity | UBAP2 Polyclonal Antibody detects endogenous levels of protein. |
| Immunogen Description | Synthesized peptide derived from part region of human protein |
| Target Name | UBAP2 KIAA1491 |
| Calculated MW | 123kD |
| Concentration | 1 mg/ml |
| Formulation | Liquid in PBS containing 50% glycerol, and 0.02% sodium azide. |
| Storage | Store at -20°C/1 year |

Application Details

WB 1:500-2000

ELISA 1:5000-20000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.