

# Rat Vascular Endothelial Growth Factor A (VEGFA) ELISA Kit

Catalog No: #EK20791

Package Size: #EK20791-1 48T #EK20791-2 96T



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Vascular Endothelial Growth Factor A (VEGFA) ELISA Kit
Applications	ELISA
Species Reactivity	Rat

## Application Details

Detect Range:15.6-1000pg/mL

Sensitivity:5.4pg/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.