

COX14 rabbit pAb

Catalog No: #62033



Package Size: #62033-1 50ul #62033-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	COX14 rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human,Mouse,Rat
Specificity	This antibody detects endogenous levels of human COX14
Immunogen Description	Synthesized peptide derived from human COX14
Target Name	COX14 C12orf62
Calculated MW	6kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.