

# Mouse Extracellular Sulfatase Sulf-1 (SULF1) ELISA Kit

Catalog No: #EK20794

Package Size: #EK20794-1 48T #EK20794-2 96T



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Extracellular Sulfatase Sulf-1 (SULF1) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse
Target Sequence	Detection Method:sandwich

## Application Details

Detect Range:78.1-5000pg/mL

Sensitivity:33.4pg/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.