

Galactose Microplate Assay Kit

Catalog No: #AS0207



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Galactose Microplate Assay Kit
Applications	Functional Assay
Species Reactivity	All
Storage	Shipped and store at 4°C for 6 months.

Application Details

Description: Detection and Quantification of Galactose Content.

Detection Method: Colorimetric

Detection Range: 0.01 mmol/L - 1 mmol/L

Sample Types: Serum, Plasma, Urine, Saliva, Milk, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids

Product Description

Detection and Quantification of Galactose Content.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.