

Monkey Thyroid Stimulating Hormone (TSH) ELISA Kit

Catalog No: #EK20800

Package Size: #EK20800-1 48T #EK20800-2 96T



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Monkey Thyroid Stimulating Hormone (TSH) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Monkey

Application Details

Detect Range:0.156-10mIU/L

Sensitivity:0.073mIU/L

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.