

Histamine Microplate Assay Kit

Catalog No: #AS0208



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Histamine Microplate Assay Kit |
| Applications | Functional Assay |
| Species Reactivity | All |
| Storage | Shipped and store at 4°C for 6 months. |

Application Details

Description: Detection and Quantification of Histamine Content.

Alternative Names: 2-(4-Imidazolyl)ethylamine; 1H-Imidazole-4-ethanamine; Ergotidine

Detection Method: Colorimetric

Detection Range: 0.01 mmol/L - 1 mmol/L

Sample Types: Tissue extracts, Other biological fluids

Product Description

Detection and Quantification of Histamine Content.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.